

Foyomed[®]

NINGBO FOYOMED MEDICAL INSTRUMENTS CO., LTD.

ADD:9F, Ningshing Building, No.666, Jiangnan Road,
Hi-Tech Zone, Ningbo, China
Tel: 86-574-81877115 Fax: 86-574-81877057
<http://www.foyomed.com> www.foyocare.com
E-mail: products@foyomed.com

Mflab[®]

NINGBO MFLAB MEDICAL INSTRUMENTS CO., LTD.

ADD: NO. 508, Yindong Road (N) , Yinzhou Economic Development zone,
315145, Ningbo, China
Tel:86-574-81877257 Fax: 86-574-81877057
<http://www.breathcare.cn> www.mflab.com
E-mail: info@mflab.com

SILICONE MEDICAL PRODUCTS

NINGBO FOYOMED MEDICAL INSTRUMENTS

Establishment

FOYOMED is established in the year 2003 as a professional import and export corporation for medical products. It is the headquarters of MFLAB Medical with production area covering 15,500 m2. Among it, it is including 100000 grade purification workshops and cleaning room, totally 9,200 m2. There are more than 250 workers in the production line.

Team Work

FOYOMED has a high efficient and professional team, full of wisdom and passion, with 15 years of working experience. We take the diligent, professional and responsible spirit, committed to provide customers good products and service.

Export Reputation

Our FOYOMED products have been exported to over 100 countries worldwide through a network of distributors and subsidiaries, covering North American, South American, European, Asian and African countries.

Principle

We FOYOMED adheres to "For Your Smile" principle to strive to achieve excellent quality, responsible price and first-rate service in operation, and cherish every business opportunity.

Leading Products

FOYOMED has a wide range of products in business with leading products respiratory products, disposable medical products, medical equipment and laboratory products.

Brands

There are five brands for our FOYOMED: Foyomed, Foyocare, FoyoLab, Mflab and Smilecare in distribution.

Quality System

FOYOMED is the headquarters of MELAB Medical which has been passed by GMP, ISO13485, US FDA on-site Audit, Saudi Arabian FDA on-site Audit, Malaysian MOH on-site Audit, TUV CE 0197 certification of more than 40 products, and a number of patents. Our FOYOMED also has been passed by ISO 13485 and CE 0123.

LB5040S Silicone Reinforced Endotracheal Tube(Cuffed)

Cat. No.	I.D.(mm)	O.D.(mm)
LB504035S	3.5	5.3
LB504040S	4.0	6.0
LB504045S	4.5	6.7
LB504050S	5.0	7.3
LB504055S	5.5	8.0
LB504060S	6.0	8.7
LB504065S	6.5	9.3
LB504070S	7.0	10.0
LB504075S	7.5	10.7
LB504080S	8.0	11.3
LB504085S	8.5	12.0
LB504090S	9.0	12.7
LB504095S	9.5	13.3



Made of medical grade silicone, better biocompatibility
High pressure silicone cuff maintains good sealing
Silicone tube with spiral reinforced wire reducing risk of kinking
100% medical grade silicone and soft tip, reduce the damage to the tracheal wall and mucosa
Single-use and reusable 2 types available, both of the two types can with stylet equipped
Reusable type is autoclave sterilization only, can be reused up to 10 times

LB3210 Disposable Silicone Laryngeal Mask

- Medical grade silicone for super compatibility and patient comfort
- Available with seven sizes for all patient use

Cat.NO.	Type	Applicable People
LB321010	1.0	< 5kg
LB321015	1.5	5-10kg
LB321020	2.0	10-20kg
LB321025	2.5	20-30kg
LB321030	3.0	30-50kg
LB321040	4.0	50-70kg
LB321050	5.0	>70kg



LB3220 Reusable Silicone Laryngeal Mask

- Reusable Silicone Laryngeal Mask are designed to be sterilized and used up to 40 times
- Durable high temperature resistant material suitable for autoclave sterilization
- Additional user card for record of device usage times and patient information is supplied in each pouch

Cat.NO.	Type	Applicable People
LB322010	1.0	< 5kg
LB322015	1.5	5-10kg
LB322020	2.0	10-20kg
LB322025	2.5	20-30kg
LB322030	3.0	30-50kg
LB322040	4.0	50-70kg
LB322050	5.0	> 70kg



LB3110
Silicone Manual Resuscitator

- ⊙ The set with oxygen reservoir bag, Silicone mask and oxygen tubing
- ⊙ 100% medical level silucone Material, Latex free
- ⊙ Additional accessories (PEEP Valve, B Valve,etc.) and customized label are available
- ⊙ All parts(excluding the oxygen tubing & oxygen reservoir)are fully autoclave up to 134°C
- ⊙ Reusable for 20 times
- ⊙ Size: Adult/Pediatric/Infant



LB3110L
Liquid Silicone Manual Resuscitator (Liquid Silicone has better elasticity)

- ⊙ The set with oxygen reservoir bag, Silicone mask and oxygen tubing
- ⊙ 100% medical level silucone Material, Latex free
- ⊙ Additional accessories (PEEP Valve, B Valve,etc.) and customized label are available
- ⊙ All parts(excluding the oxygen tubing & oxygen reservoir)are fully autoclave up to 134°C
- ⊙ Reusable for 20 times
- ⊙ Size: Adult/Pediatric/Infant





LB3030
Anesthesia Mask (one piece silicone mask)

- ⊙ All-in-one piece soft silicone cushion for better seal and comfort.
- ⊙ Ergonomic shape for lower dead space and better air circulation..
- ⊙ 100% Medical level Silicone material, Latex free.
- ⊙ Reusable all parts are fully autoclave up to 134°C
- ⊙ Available in 7 sizes.

Cat.NO.	Size	Connector (ID.)	Package
LB303000	0#,Neonate	15mm OD	Ploy bag
LB303001	1#,Infant	15mm OD	Ploy bag
LB303002	2#,Child	22mm ID	Ploy bag
LB303003	3#,Adult S	22mm ID	Ploy bag
LB303004	4#,Adult M	22mm ID	Ploy bag
LB303005	5#,Adult L	22mm ID	Ploy bag
LB303006	6#,Adult XL	22mm ID	Ploy bag



LB3031
Anesthesia Mask (rounded silicone mask)

- ⊙ Soft silicone cuff for promoting better seal and patient comfort.
- ⊙ 100% Medical level Silicone material, Latex free.
- ⊙ Rounded shape design.
- ⊙ Reusable all parts are fully autoclave up to 134°C
- ⊙ Available in 3 sizes.

Cat.NO.	Size	Connector (ID.)	Package
LB303100	0#,Neonate	15mm OD	Ploy bag
LB303101	1#,Infant	15mm OD	Ploy bag
LB303102	2#,Child	22mm ID	Ploy bag

LB3032
Anesthesia Mask (two piece silicone mask)

- ⊙ Soft silicone cushion encased in hard plastic shell for ease to handle and better seal.
- ⊙ 100% Medical level Silicone and PC material, Latex free.
- ⊙ Reusable all parts are fully autoclave up to 134°C
- ⊙ Available in 3 sizes.

Cat.NO.	Size	Connector (ID.)	Package
LB303203	3#,Adult S	22mm	Ploy bag
LB303204	4#,Adult M	22mm	Ploy bag
LB303205	5#,Adult L	22mm	Ploy bag



LB3033
Manual Sputum Processor


Cat.NO.	Size	Package
LB303301	0#	Ploy bag
LB303302	1#	Ploy bag
LB303303	2#	Ploy bag





**LB3041S
CPAP Mask**

- ⊙ Size:M,L
- ⊙ Made of high-quality liquidmedical silicone and suitable for facial features.It will not be aging, hardening but comfortable and durable
- ⊙ Sealing edge of contoue line and head belt eliminate discomfort caused by the oppression to facial nerve
- ⊙ The design of nasal part fits well witj the nose and leaves little air blind angle, which limits breath of exhaust and increases comfort
- ⊙ Comfortable, durable silicone skirt design, five points head belt, anti-fog design of standard nose and mouth cover
- ⊙ The sizes which can meet almost all type of faces and do not leak



**LB6440S
Silicone Nasogastric Tube**

- ⊙ Made of 100% medical grade silicone, the tube is soft and clear
- ⊙ X-Ray opaque line through the whole tube
- ⊙ The balloon is glued to the main tube both inside and outside, it is elastic and flexile
- ⊙ It's suitable for gastrostomy patient.

Cat. No.	Size(Fr)	Color Code	Length (mm)	Marked Positizion(cm)
LB612312	12Fr/5ml	White	240	2/4/6
LB612314	14Fr/5ml	Green	240	2/4/6
LB612316	16Fr/5ml	Orange	240	2/4/6
LB612318	18Fr/20ml	Red	240	2/4/6
LB612320	20Fr/20ml	Yellow	240	2/4/6
LB612322	22Fr/20ml	Purple	240	2/4/6
LB612324	24Fr/20ml	Dark Blue	240	2/4/6

LB6220S Silicone Suction Catheter



- ⊙ Size: Fr6-Fr20
- ⊙ Made of high-quality imported silicone rubber materials and with the method of platinum catalytic technology.
- ⊙ The tube is soft and without stimulation. There is smooth plug at the top to avoid injury caused respiratory. damage and bleeding, significantly reducing patients'pain.
- ⊙ Obvious scale marks, easy to know the depth of insertion.
- ⊙ Applicable to sputum suction in various clinical departments.




LB6410S Silicone Stomach Tube

- ⊙ Made of 100% medical grade silicone, the tube is soft and clear.
- ⊙ X-ray opaque line through the whole tube.
- ⊙ Clear depth marks on the tube.
- ⊙ Size: Fr6-Fr18

LB6420S Silicone Ryle's Stomach Tube

- ⊙ Made of 100% medical grade silicone, the tube is soft and clear.
- ⊙ X-ray opaque line through the whole tube.
- ⊙ Clear depth marks on the tube.
- ⊙ Stainless steel balls are sealed in the tip end; it is easier to indwell the tube at a certain position
- ⊙ Size: Fr6-Fr18





**LB7010S
Silicone Chest Drain**


- ⊙ Made of 100% imported medical grade silicone.
- ⊙ X-ray opaque line through the whole tube.
- ⊙ Size: Fr12-Fr36

All Silicone Foley Catheter

- ⊙ 100% silicone
- ⊙ Color-coded for easy size identification
- ⊙ Alternative types are available
- ⊙ Sterile by EO
- ⊙ For single use only
- ⊙ Hard valves only

FY0105 2-Way All Silicone Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code
FY010506	6	3	Light Red
FY010508	8	3-5	Black
FY010510	10	3-5	Grey
FY010512	12	5-10/30	White
FY010514	14	5-10/30	Green
FY010516	16	30	Orange
FY010518	18	30	Red
FY010520	20	30	Yellow
FY010522	22	30	Violet
FY010524	24	30	Blue
FY010526	26	30	Pink



All Silicone Foley Catheter

- ⊙ 100% silicone
- ⊙ Color-coded for easy size identification
- ⊙ Alternative types are available
- ⊙ Sterile by EO
- ⊙ For single use only
- ⊙ Hard valves only

FY0106 3-Way All Silicone Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code
FY010616	16	30	Orange
FY010618	18	30	Red
FY010620	20	30	Yellow
FY010622	22	30	Violet
FY010624	24	30	Blue
FY010626	26	30	Pink



Fy0122 All Silicone Foley Catheter with Temperature Sensor

- ⊙ All silicone foley catheter with temperature sensor is designed to be used before, during, and after surgery to monitor core body temperature.
- ⊙ It is made of 100% medical grade silicone, and is available with long term placement.
- ⊙ We also offer a various kind of temperature extension cable for different indication.

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code
FY012208	8	5	Black
FY012210	10	5	Grey
FY012212	12	10	White
FY012214	14	10	Green
FY012216	16	30	Orange
FY012218	18	30	Red
FY012220	20	30	Yellow
FY012222	22	30	Violet
FY012224	24	30	Blue
FY012226	26	30	Pink



Catheter (Fr)	Length (cm)	Balloon Volume (ml)
18	40	80

FY0123 Disposable Cervical Ripening Balloon

Introduce:

The Cervical Ripening Balloon is used for physical dilation of cervical canal prior while the cervix is unfavorable for induction. balloons are inflated with sterilized saline.

By using the pressure of both balloon, the cervix will synthesize and release prostaglandin, and promote the cervix ripening.

Features

1. Shorten the process of birth, reduce the pain of parturition
2. Medical grade material ensures the highest level of safety
3. Obvious effect on promoting the cervix ripening
4. Ripens and dilates the cervix without pharmaceuticals
5. Silicone body material eliminates the issue of patient sensitively to latex



Catheter (Fr)	Length (cm)	Balloon Volume (ml)
24	55	500

FY0124 Postpartum Balloon

Introduce:

Used for Temporary control or reduction of postpartum uterine bleeding when conservative management is warranted.

For while the balloon was inflated with sterilized saline, the postpartum balloon applies pressure to the uterine walls that will help to hemostasis.

Features:

1. Easy to operate and position,
2. safe and efficient way to promote hemostasis
3. Can be used for both vaginal and cesarean delivery
4. Silicone body material eliminates the issue of patient sensitively to latex



**FY0121
Silicone Condom Catheter**

- ⊙ Made of 100% imported medical grade silicone
- ⊙ Pop On: Suitable for old person or sports person.
- ⊙ Standard: Suitable for regular length and adhesive for daily wear.
- ⊙ Wide Band: More adhesive and viscous force than standard for added security.


Cat.No.	Type	Diameter				
		S	M	L	XL	XXL
FY0121P	Pop On	24mm	28mm	31mm	35mm	40mm
FY0121S	Standard	24mm	28mm	31mm	35mm	40mm
FY0121W	Wide Band	24mm	28mm	31mm	35mm	40mm



**LB6141
Mucus Trap**

Cat. NO.	Specification	Quantity (pcs)	Measurement(cm)
LB614140	Mucus trap,40cc	50	40*27*20
LB614180	Mucus trap,80cc	50	40*27*20

For sputum specimen collection, can form closed sputum specimen collection and transportation, clean, convenient, sanitary



**FY0235
Closed Wound Drainage System Kit**

- ⊙ Made of 100% imported medical grade silicone
- ⊙ This product includes 100cc/200cc Silicone Reservoir, 7mm/10mm Flat drains or other round drains and stainless steel trocar
- ⊙ Can use separately or using Trocar for insertion and combining with silicone Reservoir for fluids syction and collection
- ⊙ Used for various surgeries
- ⊙ Supplied sterile, EO sterilization




**FY0236
Silicone Reservoir/Bulb**

- ⊙ Made of 100% imported medical-grade silicone, with anti-reflyx valve.
- ⊙ Can compatible with any of following drains: Flat Drain, Fluted Drain and Round Drain.
- ⊙ Squeeze the reservoir and connect with drains for suction, the low-vacuum reservoir can make stable and effective suction for exudates.

Cat. NO.	Size(cc)	Marked Position(ml)
FY023601	100	25
FY023602	200	30
FY023603	400	40

*Clamp available.

**FY0109S
Silicone Penrose Drainage Tube**




- ⊙ 100% medical grade silicone for superior biobiocompatibility
- ⊙ X-Ray opaque line throughout the length
- ⊙ Soft and smooth tube minimizes trauma, tooth-like structure to minimize clog or obstruct
- ⊙ Clear connecting tube provides insoection and drainage visibility
- ⊙ Sterilized double package
- ⊙ Size: 6mm x 30cm, 6mm x 50cm, 8mm x 30cm, 8mm x 50cm, 10mm x 30cm, 10mm x 50cm, 12mm x 30cm, 12mm x 50cm

FY0243
Silicone Round Channel/Fluted Drain

- Made of 100% imported medical grade silicone.
- With X-ray opaque, used for closed wound suction draining.
- These intrnal four channels, or fluted, are utilized to quickly remove exudates.
- Using one-time extrusion, and hubless of the tube minimizes pain associated with device removal.

Cat. NO.	Size(Fr)	O.D.(mm)	Length(cm)
FY024310	10	3.3	120
FY024315	15	5.0	120
FY024319	19	6.3	120
FY024324	24	8.0	120

*Length could be customized.



FY0244
Silicone Round Perforated Drain

- Made of 100% imported medical grade silicone.
- X-ray opaque line through the whole tube.
- It can be equipped with reservoir and trocar.
- Clear tubing provides, inspection and drainage visibility.

Cat. NO.	Size(Fr)	O.D.(mm)
FY024410	10	3.0
FY024412	12	4.0
FY024414	14	4.7
FY024416	16	5.3
FY024418	18	6.0
FY024420	20	6.7



*Length could be customized.

**FY0245
Silicone Flat Perforated Drain**

- Made of 100% imported medical grade silicone.
- This Drain combines a transparent silicone round tube and a white flat drain.
- The inner wall of the implanted flat segment features tooth-like to prevent possible occlusion when tube is pressed or kinked.
- 3/4 or full perforation are available for optimal drainage with X-ray opaque effectiveness throughout its entire implanted segment.
- Using Trocar for insertion and combining with silicone Reservoir for fluids suction and collection.

Cat. NO.	Width(mm)	Length(cm)
FY024507	7	25
FY024510	10	30
FY024513	13	34

*Length could be customized.




**FY0242
Silicone Flat Fluted Drain**

- Made of 100% imported medical grade silicone.
- This Drain combines a transparent silicone round tube and a white flat drain. The inner wall of the implanted flat segment features multiple narrow ducts to prevent possible occlusion when the tube is pressed or kinked.
- Full fluted design is available for optimal drainage with X-ray opaque effectiveness throughout its entire implanted segment.
- Using Trocar for insertion and combining with silicone Reservoir for fluids suction and collection.

Cat. NO.	Width(mm)	Length(cm)
FY024208	8	30
FY024211	11	30


*Length could be customized.



FY0120
T-shaped Silicone Perforated Drain


- Made of 100% imported medical grade silicone.
- With X-ray Opaque and super biocompatibility, it can be used for wound draining.
- Size: Fr10-Fr20

Cat. NO.	Size(Fr)	O.D.(mm)
FY012012	12	4.0
FY012014	14	4.7
FY012016	16	5.3



LB532
Breathing Circuit for Ventilator

Made of high-quality imported silicone materials and with the method of platinum catalytic technology.
 Patent technology of integrated molding, high strength, not prone to coming off and being separated, good flexibility.
 Can endure high temperature and high pressure disinfection, reusable.
 Joints are made with injection molding with no gas leakage.
 Applicable to all kinds of imported and domestic breathing machine, anesthesia machine. Length could be customized and free options among glass, Y type joint, L-shaped connector, masks, breathing air.



LB531
Breathing Tube for Anesthesia Apparatus

Made of high-quality imported silicone materials and with the method of platinum catalytic technology.
Patent technology of integrated molding, high strength, not prone to coming off and being separated, good flexibility.
Can endure high temperature and high pressure disinfection, reusable.
Joints are made with injection molding with no gas leakage.
Applicable to all kinds of imported and domestic breathing machine, anesthesia machine. Length could be customized and free options among glass, Y type joint, L-shaped connector, masks, breathing air.

FY0310
Menstrual Cup

This menstrual cup is reusable, made of medical grade silicone, which is widely applied in medical care products such as infant nipples, medical catheters etc, because it is toxic free and allergy free, safe to use. It is designed to be folded inserted into your vagina to collect your menstrual blood. Once inside, it opens up and sits lower than tampon. It is held in place by your vagina muscles.

- Advantages:
- Keep cool and safe
 - Comfortable, clean and easy to use.
 - 100% medical grade silicone, no BPA or latex
 - Reusable, eco-friendly and economical
 - Leak-free protection for up to 10 hours at a time
 - Long-term use can reduce the risk of gynecological inflammation
 - Worry-free when travelling, swimming or exercising during menstrual period.

S: Recommended for woven who have never delivered vaginally or have delivered by caesarean section
L: Recommended for woven who have delivered vaginally

